


# Racks for ULT C75



## ORGANIZE AND OPTIMIZE THE FREEZER

Quality, efficiency and capacity at the best possible economy, is a vital objective in Nordic Lab's business philosophy. Our storage solutions include racks and boxes, all insuring an optimal handling and market low cost per sample to be stored.

Please find below specifications on the storage solution for the ULT C75 freezer model.



Storage capacity					
Box height	No. of boxes	No. of racks	Total no. of boxes	Racks without B&D* Product no.	Racks with B&D* Product no.
50 mm	8	5	40	N311210	N312210
75 mm	5	5	25	N311310	N312310
100 mm	4	5	20	N311410	N312410

\*Boxes and Dividers