

# Chest freezer ULT C300

110-120V/60Hz



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable

- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- ISO ETL, ATEX Certifications

## Specifications

### Technical

Gross volume	10,45 ft <sup>3</sup>
Net volume	10,03 ft <sup>3</sup>
Temp. range	-60°C to -86°C
Temp. probes	1/PT 1000
Compressor	Embraco
Refrigerant	R-489A
Noise level	55 dB
Power consumption 24/h	8 kWh
Power consumption	2.920 kWh / year
Voltage	110-120 V
Frequency	60 Hz

### Features

Audible & visual alarm	YES
Castors	4
Integrated lock	YES
Inner lids	3
Defrost	Manual
Battery Backup	48h
USB	NO
Data connection	NO
Digital controller with display	XR30CX
Data logging	NO
High/Low temperature alarm	YES
Porthole for external probe & CO2 backup	Ø 0,47 in
Connection for remote alarm	YES

### Material

Inner	Precoated white steel
Outer	Galvanized pre-coated steel

### Outer dimensions

Height (with castors)	34,8 in
Width	49,6 in
Depth (incl. condenser)	24,8 in
Depth (incl. condenser & handle)	27,6 in
Net weight	156 lbs

### Inner dimensions

Height	25/15,1 in
Width	43,3/35 in
Depth	17,3 in
Insulation	3,1 in
Insulation type	PUR
Internal humidity	< 5 %

## ULT C300



Product no.

N113053

## Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
2 in	6	3	18	N311200
	11	18	198	N311220
			<b>= 216</b>	
3 in	4	3	12	N311300
	8	18	144	N311320
			<b>= 156</b>	
4 in	3	3	9	N311400
	5	18	90	N311420
			<b>= 99</b>	